




Application Note

UTA 1024-MUX Adapter



Brian Pflugrath
Engineer, Verasonics Inc.

March 2018

Summary

This document describes the intended use of the Vantage UTA 1024-MUX adapter, the required configuration, programming model, assembly, and the basic sequence needed to complete a 1024 element acquisition made up of multiple transmit receive events.

Table of Contents

Summary	2
Table of Contents	2
1. Setup	3
Required Hardware.....	3
Required Software	3
2. Theory of Operation.....	4
3. Assembly and Mechanical	5
UTA best practices	5
Cover removal/replacement	5
Termination board connectors.....	6
Physical assembly images	6
4. Software	8
Apertures1024.m	8
5. Vernon Cable Termination Board	9
6. Optional Verasonics Cable Termination Board	10
Conceptual Implementation.....	10
Micro Coax Termination Pads	10

1. Setup

Required Hardware


- Vantage 256
- UTA 1024-MUX Adapter

Required Software

- Vantage Release 3.2.2 or later

2. Theory of Operation

The 1024-MUX adapter supports 4-to-1 multiplexing of 256 system channels to provide connection up to 1024 elements. The user can define custom apertures using the aperture helper function within the limitations of mux exclusivity. For example, system channel 1 can be connected to elements 1, 257, 513, or 769. The system channel cannot connect to more than one mux position at a tie. See diagram below.




3. Assembly and Mechanical

UTA best practices


- * 1024-MUX UTA assemblies should be handled with Electro Static Discharge (ESD) safe practices to avoid electrical component damage
- * Transducers should ONLY be assembled when the 1024-MUX UTA is unplugged from the Vantage system
- * Extreme caution should be used when plugging and unplugging termination boards so as not to damage the soldered micro coax connections
- * Plugging and unplugging of termination should be minimized to preserve connector integrity/longevity

Cover removal/replacement


To remove or replace lid locate the 8 M3x8mm screws around the perimeter of the metal lid



Termination board connectors



Physical assembly images






4. Software

computeUTAMux1024.m

The computeUTAMux1024.m helper function requires an aperture or an array of apertures to populate the appropriate channel-mapping and mux programming structures. To create one or more apertures a user must create a table of ones and zeros 1024 columns wide. computeUTAMux1024.m will check the input for correctness and will error out and notify the user if incorrect. If the transducer does not directly map from element to channel, the user must define a remap table for the adapter/system to use. Additional details can be found in the computeUTAMux1024.m function in the Utilities folder.


5. Vernon Cable Termination Board

Vernon 1024 element transducers are supplied with 8 termination boards compatible with the Vantage1024-MUX adapter. Each board connects to 128 elements. The default Vernon element to channel mapping is a direct map.




6. Optional Verasonics Cable Termination Board

Conceptual Implementation



Micro Coax Termination Pads



Pad order from left to right is 1 to 52. There are five sets of landing pads from top of the board towards the bottom which are listed as MCLP1 thru 5 below (Multi-Coax-Land-Pattern).

MCLP1

1	EL_13	14	EL_65	27	EL_141	40	EL_193
2	EL_12	15	EL_66	28	EL_140	41	EL_194
3	EL_11	16	EL_67	29	EL_139	42	EL_195
4	EL_10	17	EL_68	30	EL_138	43	EL_196
5	EL_9	18	EL_69	31	EL_137	44	EL_197
6	EL_8	19	EL_70	32	EL_136	45	EL_198
7	EL_7	20	EL_71	33	EL_135	46	EL_199
8	EL_6	21	EL_72	34	EL_134	47	EL_200
9	EL_5	22	EL_73	35	EL_133	48	EL_201
10	EL_4	23	EL_74	36	EL_132	49	EL_202
11	EL_3	24	EL_75	37	EL_131	50	EL_203
12	EL_2	25	EL_76	38	EL_130	51	EL_204
13	EL_1	26	EL_77	39	EL_129	52	EL_205

MCLP2

1	EL_26	14	EL_78	27	EL_154	40	EL_206
2	EL_25	15	EL_79	28	EL_153	41	EL_207
3	EL_24	16	EL_80	29	EL_152	42	EL_208
4	EL_23	17	EL_81	30	EL_151	43	EL_209
5	EL_22	18	EL_82	31	EL_150	44	EL_210
6	EL_21	19	EL_83	32	EL_149	45	EL_211
7	EL_20	20	EL_84	33	EL_148	46	EL_212
8	EL_19	21	EL_85	34	EL_147	47	EL_213
9	EL_18	22	EL_86	35	EL_146	48	EL_214
10	EL_17	23	EL_87	36	EL_145	49	EL_215
11	EL_16	24	EL_88	37	EL_144	50	EL_216
12	EL_15	25	EL_89	38	EL_143	51	EL_217
13	EL_14	26	EL_90	39	EL_142	52	EL_218

MCLP3

1	EL_39	14	EL_91	27	EL_167	40	EL_219
2	EL_38	15	EL_92	28	EL_166	41	EL_220
3	EL_37	16	EL_93	29	EL_165	42	EL_221
4	EL_36	17	EL_94	30	EL_164	43	EL_222
5	EL_35	18	EL_95	31	EL_163	44	EL_223
6	EL_34	19	EL_96	32	EL_162	45	EL_224
7	EL_33	20	EL_97	33	EL_161	46	EL_225
8	EL_32	21	EL_98	34	EL_160	47	EL_226
9	EL_31	22	EL_99	35	EL_159	48	EL_227
10	EL_30	23	EL_100	36	EL_158	49	EL_228
11	EL_29	24	EL_101	37	EL_157	50	EL_229
12	EL_28	25	EL_102	38	EL_156	51	EL_230
13	EL_27	26	EL_103	39	EL_155	52	EL_231

MCLP4

1	EL_52
2	EL_51
3	EL_50
4	EL_49
5	EL_48
6	EL_47
7	EL_46
8	EL_45
9	EL_44
10	EL_43
11	EL_42
12	EL_41
13	EL_40

14	EL_104
15	EL_105
16	EL_106
17	EL_107
18	EL_108
19	EL_109
20	EL_110
21	EL_111
22	EL_112
23	EL_113
24	EL_114
25	EL_115
26	EL_116

27	EL_180
28	EL_179
29	EL_178
30	EL_177
31	EL_176
32	EL_175
33	EL_174
34	EL_173
35	EL_172
36	EL_171
37	EL_170
38	EL_169
39	EL_168

40	EL_232
41	EL_233
42	EL_234
43	EL_235
44	EL_236
45	EL_237
46	EL_238
47	EL_239
48	EL_240
49	EL_241
50	EL_242
51	EL_243
52	EL_244

MCLP5

1	N/C
2	N/C
3	EL_64
4	EL_63
5	EL_62
6	EL_61
7	EL_60
8	EL_59
9	EL_58
10	EL_57
11	EL_56
12	EL_55
13	EL_54

14	EL_53
15	EL_117
16	EL_118
17	EL_119
18	EL_120
19	EL_121
20	EL_122
21	EL_123
22	EL_124
23	EL_125
24	EL_126
25	EL_127
26	EL_128

27	EL_192
28	EL_191
29	EL_190
30	EL_189
31	EL_188
32	EL_187
33	EL_186
34	EL_185
35	EL_184
36	EL_183
37	EL_182
38	EL_181
39	EL_245

40	EL_246
41	EL_247
42	EL_248
43	EL_249
44	EL_250
45	EL_251
46	EL_252
47	EL_253
48	EL_254
49	EL_255
50	EL_256
51	N/C
52	N/C

This document is intended for use by Verasonics customers only.

Do not duplicate or distribute without permission.

Verasonics, Inc.

11335 NE 122nd Way, Suite 100

Kirkland, WA 98034

USA